EL-WiFi-ULT+

High Accuracy, Ultra Low Temperature Cryogenic Vaccine WiFi Data Logger

- Measurement range -200 to +100°C (-328 to +212°F)
- High accuracy T-type probe on 2m cable
- Wirelessly stream and view data on the EasyLog Cloud, App or on a PC
- View and analyse multiple sensors, including graphing of historic data
- Configurable high and low alarms with indicator
- Easy sensor set-up using free PC software or App

The EL-WiFi-ULT+ is designed to monitor vaccines in cryogenic dry ice storage.

Data is streamed wirelessly over any compatible WiFi network and can be viewed on a PC using free software or on the EasyLog Cloud or App.

To configure the sensor for use on a given wireless network, either connect it via USB to a PC running EasyLog WiFi software on the network, or configure wirelessly using the EasyLog Cloud app on a mobile phone with access to the network. The sensor can then be placed anywhere within range of the network. If the sensor temporarily loses connectivity with the network, it will log readings until it is able to communicate again with the PC application or EasyLog Cloud (max 30 days at 10 second sample interval).

The sensor is IEEE 802.11bgn* (2.4GHz) compliant, supports WEP, WPA/WPA2 encryption and enterprise networks (PEAP, TTLS, FAST).

EL-WiFi-ULT+ has a protection rating of IP40. The unit is freestanding, but can be attached to a wall or surface using the bracket provided. The unit can be clipped in and out of the bracket as required.

SPECIFICATIONS

	Minimum	Typical	Maximum	Unit	
Battery life		>6		Months	
Temperature measurement range	-200 (-328)		+100 (+212)	°C (°F)	
USB supply voltage (@500mA)	4.5	5.0	5.5	Vdc	
Operating temperature range	-20 (-4)		+60 (+140)	°C (°F)	
Logging period (user configurable)	10 sec	10 min	12 hrs		
Transmission period (user configurable)	1 min	1 hr	24 hrs		
Temperature measurement resolution		0.1 (0.2)		°C (°F)	
Temperature display resolution		0.1			
Accuracy		±1.8 (±3.2)		°C (°F)	
IP Rating	IP40 (Bung fully inserted, not permanently powered, thermo- couple connector fitted and fully inserted, device mounted vertically.)				
Dimensions	82 x 70 x 23mm (3.22 x 2.75 x 0.91")**				
Probe dimensions	300 x Φ 1.5mm, 2m cable (11.8 x Φ 0.05, 39.3" cable)				

ACCESSORIES

PSU-5VDC-USB-USA	USB Mains Power Adapter for USA			
PSU USB-UK	USB Mains Power Adapter for UK			
PSU USB-EU	USB Mains Power Adapter for EU			
EL-P-TC-T-ULT+	Replacement T-type Thermocouple probe			
INCLUDED IN THE BOX				
EL-WIFI WALL BRACKET	Wall mounting bracket for EL-WiFi sensors			
EL-P-TC-T-ULT+	Thermocouple probe			



* MAC Address starting 98:8B:AD:2...... only

** Excluding probe and mounting bracket





CALIBRATION CERTIFICATES NOW AVAILABLE

Lascar now offers a Traceable Calibration Certificate Service on Temperature Data Loggers. Using reference equipment which has been calibrated by a UKAS/NIST accredited laboratory and using apparatus traceable to national or international standards. For more information please see www.lascarelectronics.com.



EL-WiFi-ULT+

High Accuracy, Ultra Low Temperature Cryogenic Vaccine WiFi Data Logger

EasyLog Software

Lascar's WiFi software* is available as a free download from: www.lascarelectronics.com/software/easylog-wifi. Easy to install and use, EL-WiFi-WIN allows easy connection of sensors to a WiFi network. The user can select where data is stored - the PC or the Cloud.

When setting up your logger, make sure you select the T-Type probe.

EasyLog Cloud Your Data. Anytime. Anywhere.

EasyLog Cloud harnesses the power of IoT to automate data logging and alert notifications, enabling you to monitor and manage multiple data logging devices in different locations completely remotely. The system easily scales to meet your needs. Perfect for compact systems with just a few measuring points, or corporate solutions with thousands of devices around the globe.

You will need to create an account at www.easylogcloud.com before setting up vour cloud-connected data logger.

Features at a glance*



Store your data logging records securely on the Cloud



Connect Connect data multiple users loggers from with variable multiple sites account in a single privileges account



Easily access your most important data, anywhere



Remotely manage all of your data logging devices



Never miss a critical event with flexible advanced notifications

Review and analyse your data with powerful graphing functionality

Keep track of data events and system activity with a detailed event log

.....

*Features depend on account type.

Battery Life and Power Supply

The battery can be recharged (unit must be between 0 - 40°C) via a PC, a USB +5V wall adapter, or a portable USB battery pack using the USB lead provided. It can also be permanently powered by a USB wall adapter or USB battery pack. Readings may be affected while the internal battery is being charged. However, once charged, continued connection of the charger will have no effect.

Battery life is dependent on: transmission period, WiFi encryption method, WiFi encryption key rotation frequency (determined by the router/access point), signal strength between router/access point and WiFi device, presence volume and type of WiFi traffic from other devices, sample rate and operating temperature.

Specifications liable to change without prior warning

*Requires Windows 7, 8.1, 10











EasyLog

electronics Issue 2 06/2024 Page 2 of 2